











Upper arm models										
Model No.	UA-1200BLE	UA-1020 /1020-W	UA-767F	UA-767S /767S-W	UA-704	UA-651 BLE	UA-651Plus /651SL Plus	UA-611Plus	UA-651 /651SL	UA-611
										
Automatic inflation	•	•	•	•		•	•	•	•	•
One-touch measurement	•	•	•	•		•	•	•	•	•
3-line display		•	•	•		•	•	•	•	•
Memory storage	100 x 5	90	60 x 4	60		30	60	60	30	30
Average reading of stored data	• (App)	•	•	•		• (App)	•	•	•	•
IHB indicator	IHB	IHB	%IHB	%IHB	IHB	IHB	IHB	IHB	IHB	IHB
WHO Blood Pressure classification indicator	• (App)	•	•	•		• (App)	•	•	•	•
AC inlet	USB	•	•	•		•	•		•	
Clock and calendar	• (App)	•	•			• (App)				
AM/PM memory recall w/average	• (App)	•				• (App)				
TriCheck™		•								
MyPressure		•								
Cuff fit error indicator	•	•	•	•						
Movement error indicator	•	•	•	•						
Data communication	Bluetooth® Low Energy					Bluetooth Low Energy				
Optimized inflation	•*	•	•	•		•	•	•	•	•
Optimized deflation		•					•	•		
Cuff (Standard accessory)	Wide range 22 cm - 42 cm	SmoothFit 17 cm -32 cm / Wide range 22 cm - 42 cm	SlimFit 22 cm - 32 cm	SlimFit 22 cm - 32 cm / Wide range 22 cm - 42 cm	SlimFit 22 cm - 32 cm	SlimFit 22 cm - 32 cm	SlimFit (ISO) 22 cm - 32 cm /23 cm - 37 cm	SlimFit (ISO) 22 cm - 32 cm	SlimFit 22 cm - 32 cm /23 cm - 37 cm	SlimFit 22 cm - 32 cm
Clinical validation paper	ISO81060-2	BHS A/A	ESH Pass	ESH Pass	BHS A/A	ESH Pass	BHS A/A	BHS A/A	ESH Pass	ESH Pass
Clinical test report ISO 81060-2:2013	•	•	•	•	•	•	•	•	•	•

*Measurement during inflation