







WIRELESS UPPER ARM BLOOD PRESSURE MONITOR



KEY FEATURES:

-  Wireless
-  One button operation
-  Irregular Heartbeat Detection
-  Pressure Rating Indicator on app
-  Unlimited memory on app
-  5 year warranty



CLINICALLY
VALIDATED

A wireless blood pressure monitor that automatically records and tracks measurements via mobile app to enable trending and sharing.

Trusted Accuracy

- Clinically validated for accuracy
- Irregular Heartbeat (IHB) alerts the user to irregular heartbeat patterns
- Pressure Rating Indicator classifies blood pressure in accordance with the medically approved standards, available on app

User Friendly

- Easy, one button operation for quick and accurate measurement
- Send your measurements wirelessly via Bluetooth to your supported phone/tablet
- Free companion A&D Heart Track app automatically stores readings, tracks progress with simple graphs and allows you to share results with others
- 30 reading memory in off-line mode, unlimited on app

- Easy to apply medium cuff that comfortably fits most arms (9-14.6")
- Gentle inflation enables a more comfortable measurement

The A&D Promise

- Professional grade performance you can trust
- Simple integration to telehealth systems with full service technical support
- 5 year warranty on monitor, 2 year warranty on cuff
- Cuffs are not made with natural rubber latex
- Battery life indicator



Technical Specifications

Measurement Method	Oscillometric
Memory	Unlimited on App
Cuff Size	9-14.6" / 23-37 cm
Measurement Range	Pressure: 0-299 mmHg Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg Pulse: 40-180 beats / minute
Measurement Accuracy	Pressure: ±3 mmHg Pulse: ±5%
Power Supply	4 x 1.5V batteries (R6P, LR6 or AA) or AC adapter (TB-233)
Number of Measurements	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries)
Monitor Dimensions	3.8 x 2.7 x 5.1 in / 9.6 x 6.8 x 13 cm
Monitor Weight	8.8 oz without batteries
Display Size	1.4 x 2.1"
Clinical Test	ISO81060-2: 2013
Operating Conditions	+50°F to +104°F (10°C to 40°C) / 15% to 85% RH / 800 to 1060 hPa
Storage/Transport Conditions	-4°F to +140°F (-20°C to 60°C) / 10% to 95% RH / 700 to 1060 hPa

Packaging and Shipping Specifications

	Dimensions (WxDxH) (in)	Weight (lb)	Individual Item Cube (cubic feet)	UPC / GTIN
UA-651BLE Individual Unit (Box)	5.12 x 7.64 x 4	1.52	0.09	0-93764-60332-7
UA-651BLE Case (10 units)	21.65 x 11 x 8.46	15.87	1.18	1-00-93764-60332-4

What is In The Box

- 1 Upper arm blood pressure monitor
- 1 Blood pressure cuff (9-14.6")
- 1 Instruction manual
- 1 Quick start guide
- 1 Carrying case
- 4 AA batteries

