

CONNECTED MONITORING

OUR FOCUS ON THE FUTURE OF TECHNOLOGY IS ENGRAINED IN THE FOUNDING PRINCIPLES OF A&D

A&D is a pioneer in connected health, developing the first connected monitor and weight scale almost 15 years ago for clinical application use.

Our pursuit of advanced and intuitive technology empowers consumers to track their health data accurately and manage their health simply.



AHA/ACC NEW GUIDELINES

Hypertension classified as 130/80 will categorize 46% of the population as hypertensive. The change to guidelines highly impacts younger patients, women with hypertension under 45 will double while men with hypertension under 45 will triple.

Lifestyle modifications such as diet exercise and **home blood pressure monitoring** are recommended as the first line of treatment.

$$\frac{140}{90} \rightarrow \frac{130}{80} \begin{matrix} \text{(sys)} \\ \text{(dia)} \end{matrix}$$



SIMPLE CONTROL

One Button Operation

Simply push one button to take a measurement. The reading will be displayed on the screen in large and easy-to-read numbers.


Easy Sync

By downloading the **free** app, patients will no longer need to keep a cumbersome logbook. The connected blood pressure monitors store readings so patients have the option of periodically syncing to the app.

A CONNECTION

To Doctors and Healthcare Practitioners

Connectivity is designed to improve patient compliance. Doctors have the ability to adjust patient care based on supportive data and identify trends in the app interface without having to rely on inaccurately kept or messy logbooks.



To Family and Caretakers

We understand that some patients are lucky enough to have a supportive network, the app is designed with the caretaker in mind. As a patient advocate, they need the right data to impact important medical decisions.

CONNECTED HEALTH SOLUTIONS

A&D provides a range of connected products providing a window to your health through our app.

TRACK YOUR

Blood Pressure SpO2 Glucose Weight

HOW TO PAIR

1



Download the free "A&D Connect" App. Open the App and turn device monitor on

2




Select device on phone. Press and hold 'START' on monitor

3



Once 'Pr' is displayed on monitor, select 'pair' on mobile device

4



When 'END' is displayed on monitor, pairing is complete

	Blood Pressure Category	Systolic (mmHg)		Diastolic (mmHg)
	NORMAL	<120	&	<80
	ELEVATED	120 – 129	&	<80
	HYPERTENSION STAGE 1	130 – 139	or	80 – 89
	HYPERTENSION STAGE 2	140+	or	90+
	HYPERTENSIVE CRISIS	180+	&/or	120+